
2-in-1 Digital Microscope

User Manual V1.0


Note: This device only supports iOS (iPhone) devices and Android devices with OTG functionality. It does not support any computer systems!

Contents

1.Read Me First	03
2.Product Introduction	04
3.Product Description	05
4.How to Use	06
4.1.App Download	06
4.2.How to Connect with Cell Phone	07
5.APP Introduction	08
6.Product Parameters	09
7.FAQ	10

1.Read Me First

① Please open the lens cap before use.


② This device combines digital microscope and optical microscope. The specific magnification effect depends on the actual captured images.


③ The optimal focal distance for the microscope is 10-40mm, different object distances correspond to different magnifications. You need to adjust the focus wheel to achieve a clear view of the object.

2.Product Introduction


This device is a high-definition smartphone direct-connect microscope. When you use the microscope, it can be directly connected to the smartphone's interface, allowing image display on the App. It supports take photos and videos, and is compatible with both iOS and Android devices with OTG functionality. It can easily be used in various fields, including:

Textile Industry

Print Inspection

PCB Precision Inspection


Educational Purposes


Hair Inspection

Skin Inspection

Biological Observation


Jewelry and Coin
(Collectibles) Inspection

Visual Assistance

3.Product Description

1 Lightning (iPhone)

2 Type-C

3 Light Adjustment Button

4 Photo Button

5 Focus Wheel

6 Led Lights

7 Lens Cap

4.How to Use

Compatible with iOS devices and Android devices with OTG functionality.

4.1 App Download

Method 1

For iOS users: Please search for "**UseePlus**" at App Store for download and install it.

For Android users: Enter the following URL in your mobile browser to download:


<https://active.clewm.net/FpZAiR?>

Method 2

Scan the QR code, click download and install:

QR code


The iOS / Android APP Icon

4.2 How to Connect with Cell Phone

Based on your phone type, choose the corresponding interface to connect the device to your phone and you will see the device's LED light turn on.

(Note: If the device's light doesn't turn on, please check the interface is properly connected. If you're using an Android device, please go to the settings page to verify if your device's OTG function is enabled.)

Type-C or Lightning (for iPhone) interface

LED lights

5.APP Introduction

This screenshot shows the home screen of the app. At the top, there is a status bar with the time 1:43, signal strength, and 5G connectivity. Below the status bar, there is a large circular graphic with a checkmark and the text '设备已连接' (Device connected). At the bottom, there are two main menu items: '摄像头' (Camera) and '相册' (Album). A settings icon is located in the top right corner.

1 Settings

2 Connection status prompt

3 Enter the imaging interface

4 Enter the album

This screenshot shows the camera interface. At the top, there is a status bar with the time 1:44, signal strength, and 5G connectivity. Below the status bar, there is a back arrow on the left and three icons on the right: a rotate screen icon, a color and black & white mode icon, and an image contrast icon. The main area shows a blue textured background. At the bottom, there are three icons: a magnification icon, a photo/video icon, and a magnification status icon.

5 Back

6 Rotate the screen

7 Color and black & white mode

8 Image contrast

9 Magnification

10 Photo/Video

11 Enter the album

This screenshot shows the album interface. At the top, there is a status bar with the time 2:09, signal strength, and 5G connectivity. Below the status bar, there is a back arrow on the left, a '相册' (Album) title, and a PDF icon on the right. The main area shows a grid of photos. At the bottom, there is a text prompt '已经全部加载完毕' (All loaded). A video playback icon is located in the top right corner.

12 Create PDF documents

13 Video playback

14 Photos playback

6.Product Parameters

Item	Specification
Product Name	X6D
Lens Optic Size	1/6.5"
Photo Resolution	1920*1440
Video Resolution	1920*1440
Video Format	MP4
Picture Format	JPEG
Focus Mode	Manual
Magnification	500/1000/1600
Digital Zoom	Supported
Illuminate	8 LED Lights （brightness adjustment）
Focus Range	10 to 40 mm （Viewable distant objects）
PC Connection	Not Supported
White Balance	Auto
Exposure	Auto
Operating System	Android 6.0 and above, iOS 10.0 and above
Languages	7 languages
Lens Structure	2G+IR
Aperture	F4.5
Lens Angle	16°
Interface and Signal Transmission Mode	Type C, Lightning
Operating Temperature	-10°C ~ +55°C
Operating Humidity	30%~85%Rh
Operating Current	110mA
Power Dissipation	0.6W

7.FAQ

Q1: Which systems is this device compatible with?

A1: Android 6.0+ (with OTG functionality) / iOS 10.0 and above.

Q2: Why is the app not displaying images?

A2.1: Check if the LED light is illuminated. If it's not, it may be due to a poor interface connection. Please unplug and reinsert the device.

A2.2: For Android users, please check if the OTG function is enabled in your device settings.

Q3: Why is the image very blurry when viewing?

A3.1: Please check if you have adjusted the microscope to the optimal viewing distance. You can achieve the clearest focus by slowly rotating the focus wheel left and right.

A3.2: Please check if you have opened the lens cap in front of the lens.

Q4: What is the reason for the static black spots displayed on the screen when using the device for viewing?

A4: The static black spots may be due to dust entering the device. You can try tapping the product gently to see if it dislodges the dust.


Q5: What should I do if the product's imaging is lagging?

A5: When the product's imaging is lagging, please make sure your phone has at least 20% battery life remaining before using the device.

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