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

**Model: CM4**

## **User Manual V1.0**

Please read this instruction manual carefully  
and keep it for future reference.


# Contents

<b>1. Parts and Functions</b>	<b>03</b>
• Parts and functions	03
<b>2. Interface Introduction and Usage</b>	<b>05</b>
<b>A.Mode introduction</b>	<b>05</b>
• Photo mode	05
• Video mode	05
• Playback mode	06
<b>B.Menu interface</b>	<b>07</b>
• Menu mode	07
• Function description	08
<b>3. Work with PC</b>	<b>09</b>
• Software download	09
• Operation steps	10
• USB mode introduction	11
<b>4. Product Parameters</b>	<b>12</b>
<b>5. Notice</b>	<b>13</b>

---

# 1. Parts and Functions

The plan used in this manual is for reference only, please refer to the actual plan.

## • Parts and functions

NO.	Function
①	Left key
②	Reset key
③	Menu key
④	Right key
⑤	Type-C interface
⑥	TF card slot
⑦	LED light
⑧	Indicator light
⑨	Display screen
⑩	Power key
⑪	Focus wheel

**Left key:** Enter the Menu or Playback interface , you can select the previous item ; When entering the Photo and Video mode , you can zoom in / out.

**Right key:** Enter the Menu or Playback interface , you can select the next item ; When entering the Photo and Video mode , adjust the brightness of the light .

**Menu key:** Short press to switch three modes: Photo / Video / Playback ;Long press to enter the Menu interface ;short press again to exit the Menu or return to the previous level .

**Indicator light:** Red indicator light is on when charging ;Green indicator light is on after the full charge ;Blue indicator light is always on when starting up .

**Power key:** Long press for 3 seconds to turn on / off the microscope; Short press to confirm .

**Display screen:** 2.0-inch IPS display screen .

**Focus wheel:** Rotate this wheel to focus the observation object and slowly adjust to a clear picture .

**TF card slot:** Support up to 32G memory card.(Note: Please choose the C10 high-speed card .)

**Type-C interface:** Standard USB 2.0 interface, which supports charging and data transmission .

**Reset key:** When the microscope is abnormal, intercept the hole with a fine tip object,which can be forced to shutdown .

**LED light:** when the ambient light is insufficient, the light can be replenished .


## 2.Interface Introduction and Usage

[ ▲ Note:This product can save 58 photos in its own memory without inserting a TF card. Photos of its own memory cannot be read by other devices for security reasons.]

### A. Mode introduction


#### ● Photo mode

When the icon in the upper left corner of the microscope is "📷", it is currently in the Photo mode; in the boot state, short press the Menu key on the microscope to switch to the Photo mode.


#### ● Video mode

When the icon in the upper left corner of the microscope is "📹", it is currently in the Video mode; in the boot state, short press the Menu key on the microscope to switch to the Video mode.

## • Playback mode

When the icon in the upper left corner of the microscope is "▶", it is currently in the Playback mode; in the boot state, short press the Menu key on the microscope to switch to the Playback mode.


As shown, enter the Playback mode to view the pictures and videos,  
When viewing the video, press Power key to play / pause;  
When viewing the current picture / video, press Left / Right key to  
turn the page, and press Menu key again to return to Photo mode.


## B. Menu interface

### • Menu mode

Long press Menu key on the microscope and enter the Menu interface in Photo / Video mode as shown in the following figure:

Enter the Menu when viewing the file in Playback mode, as shown in the figure below (you can Delete Current / Delete All the picture or video).

## ● Function description

**Image Resolution :** 2M, 3M, 5M, 8M, 10M, 12M, 16M.

**Video Resolution :** 1080P, 720P, VGA

**Date Stamp :** When opened, the date and time mark will be displayed in the lower right corner of the photo/video.

**Language :** Set up languages (English, Russian, French, Dutch, Traditional Chinese, Simplified Chinese, Japanese, Korean, German, Portuguese, Polish, Arabic, Italian, Spanish).

**Brightness:** There are 8 brightness options.

**Auto Power Off :** After selecting the automatic shutdown time, the microscope will automatically shut down when unmanned; no automatic shutdown is required, please select "Off".

**Screen Savers :** Select the appropriate time to make the device automatically turn off the screen when no one operation.

**Format :** This action removes all data in the TF card.  
( ▲Note: It cannot be restored after deletion, Please caution).

**Date/Time :** The date / time parameter displayed by the microscope can be set.

**Default Setting :** All settings will return to the factory default parameters.

**Version :** View the current system firmware version.

## 3.Work with PC

### • Software download

Note : The following three methods can be downloaded.Please select the applicable version according to the computer system. If the measurement function is not required, the own Camera software of the computer system can also be used normally.

①Please download the applicable software from: <http://www.hvscam.com/soft.asp> (for: Windows 7, Windows 8.1, Windows 10, Windows 11 and Mac).

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Support: W01 DM1 DM2 DM3 DM4 DM7 DM8 DM9 MS1 MS2 MS3 MS4 MS5

 Click to read the instructions


② For Windows 8 or higher system, search "hiView" () in the Microsoft Store to download it and you can also be viewed in Camera() .

③ For Mac OS X 12 or higher system, search "hiView" in the APP Store to download it and you can also be viewed in Photo Booth() .

## • Operation steps


### Step 1 :

Connect to the computer and the microscope via the Type-C data cable;


### Step 2:

Turn on the microscope, the pop-up window appears, select the "PC Camera" mode;


### Step 3:

Open the Hiview software, you can view the image, the photo/video file is automatically saved to the computer disk. (Note: When connecting the computer to view the image, please confirm that you are using a cable with data transmission function).

## ● USB mode introduction

When the device is connected to the computer via a Type-C cable, the device has three Memory / PC Camera / Charging mode for users to choose from.

### ① **Memory :**

Entering the microscope memory from the computer side, you can view and process the pictures and videos in the current TF card. At this time, the microscope only Power key and Right key are available.

### ② **PC Camera:**

After selecting this mode, a larger multiple image can be obtained through the computer camera software. At this time, the microscope Power key, Right key and Focal wheel are available.

### ③ **Charging mode:**

When the power of the microscope is low, the microscope can be used while charging. The microscope can be connected to the computer terminal with the Type-c data cable. When the microscope is turned on, select the charging mode, you can enter the charging state.

(Note: You can also use a daily mobile phone charger, computer and power bank to charge the microscope.)

## 4.Product Parameters

Item	Parameter
Product name	CM4
Maximum video pixel	1080P
CMOS size	1/6.5'
Picture resolution	2M,3M,5M,8M,10M,12M,16M
Video resolution	1080P,720P, VGA
Picture format	JPG
Video format	AVI
Magnification	500/1000
Light	8 LED lights (adjustable brightness)
Focus range	10-40mm (viewable vision)
PC resolution and image transmission rate	640*480/25 fps
PC operating system	Windows xp,Win7, Win8.1, Win10 Mac OS x 12 and above
Languages	14 languages
Interface and signal transmission mode	USB 2.0
Operating current	250mA
Power	1.25W

## 5.Notice

- ◆ The TF card should be selected with a transmission speed of C10. If the card with a transmission speed lower than C10 is used, there may be faults such as photo stuck or damaged video files.
- ◆ Do not place the microscope in direct sunlight and avoid a rainy or wet environment to avoid equipment failure.
- ◆ Do not place the microscope near the heat source or fire source, such as electric furnace, microwave oven, oven, furnace and other environments that may produce high temperature.
- ◆ Do not modify the battery, insert foreign objects, immerse in water or other objects to avoid overheating, fire or explosion.
- ◆ Please avoid storing the microscope in a place with the risk of falling, which may cause the microscope internal lens deviation or damage to other components, and other irreparable damage.