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Model: DM4D

User Manual V1.0


Please read this user manual carefully and keep it for future reference.

Contents

1. Parts and Functions	03
1.1-Functions introduction	04
2. Interface Introduction	05
2.1-Photo mode	05
2.2-Video mode	05
2.3-Playback mode	06
2.4-Menu mode	07
3.Function Introduction	08
4. Work with PC	09
4.1-Software download	09
4.2-Operating steps	10
4.3-USB function introduction	11
5. Specifications	12
6. Warning	13

1. Parts and Functions

* **Notice:** The picture is for reference only, please refer to the actual object

1. Parts and Functions

1.1 Functions introduction

Confirm/OK button: When Photo mode, short press to take pictures and store automatically, When Video mode, short press to start/end recording; When viewing video files in Playback mode, short press to play video; When entering the Menu interface, short press to confirm the selection.

Left/Right button: When entering the Menu or Playback mode, it's the selection button; When entering the Photo or Video mode, it's the zoom in/out button.

Menu/Switch button: Long press to enter the Menu interface, and short press again to exit the Menu or return to the previous level; Short press to switch between three modes of Photo/Video/Playback.

Indicator light: When charging, the red indicator light is on; When fully charged, the green indicator light is on; When powering on, the blue indicator light is on.

Power button: Long press for 3s to turn the microscope on/off.

TF card slot: Support maximum 32G memory card (Note: please use C10 high-speed card).

Type-C interface: This interface supports charging and data transfer.

Focus wheel: Rotate this wheel to focus on the observation object, slowly adjust until the picture is clear.

Reset button: When the microscope is abnormal, poke this hole with a fine pointed object to force the power off.

Light wheel: Turn the wheel to adjust the LED light brightness.


LED lights: When the observation environment is not enough light, it can make up light.

2.Interface Introduction

[Note: Before using the microscope, to ensure better storage of the photos /videos, please insert the TF card into the microscope first.]

2.1-Photo mode


When the icon in the upper left corner of the microscope is "📷", it is in Photo mode, When the microscope is on, short press the Menu button on the microscope, it can switch to Photo mode.

2.2-Video mode


When the icon on the upper left corner of the microscope is "📹" it is in Video mode, When the microscope is on, short press the Menu button on the microscope, it can switch to Video mode.

2.Interface Introduction

2.3-Playback mode

When the icon in the upper left corner of the microscope is "" it is in Playback mode, in the power on state, short press the Menu button on the microscope, it can switch to Playback mode.

As shown in the figure, enter the Playback mode to view pictures and video files, view the pictures selected by pressing the OK key to view, View the video by pressing the OK key for play/pause; View the current picture/video file, you can press the Left/Right selection key to turn the page, or short press the Menu button again, you can return to the previous level to re-select the files you need to view.


2.Interface Introduction

2.4-Menu mode

Long press the "M" button on the microscope to enter the menu interface in the Photo/Video mode, as shown in the following figure:

Enter the menu while viewing a file in Playback mode, as shown in the figure (you can delete/lock the current file)

3.Function Introduction

Image resolution: 1920×1080, 1280×720, 640×480

Video Resolution:1080P, 720P, VGA

Language setting:Setting language (English,Traditional Chinese,Simplified Chinese,Japanese,French,German,Korean,Italian,Portuguese,Russian,Spanish,Polish).

Time Setting:You can set the date/time parameter of the microscope.

Screen savers:Select the appropriate time to turn off the screen automatically when the unit is not operated.

Frequency: Utility frequency setting.

Auto Power-off:Select the appropriate time to make the device shutdown automatically.

Format:This operation will delete all the data in the TF card. (Note: after deletion will not be able to recover, please operate with caution).

Reset System:Restore factory settings.


Version: Check the current firmware version.

4. Work with PC

4.1-Software download

*Note: The following 3 methods can be downloaded, please choose the appropriate version according to your computer system.

- ①Please download the applicable software from the following website: <http://soft.hvscam.com>(Suitable for: Windows 7, Windows 8.1, Windows 10, Windows 11 and Mac)

- ②When the computer system is Windows 8 or higher, search for "HiView" in the "Microsoft Store" to download (as shown in the picture" ); you can also use the built-in "Camera" of the computer to observe (as shown in the picture" ).
- ③For Mac OS x 12 or above,search for “HiView” in “App Store” you can also use the built-in "Photo Booth" of Mac to observe (as shown in the picture" ).


4.Work with PC

4.2-Operating steps

Step 1: Connect the computer and microscope through the Type-C data cable;


Step 2: When the microscope is turned on and a pop-up window appears, select "PC Camera" mode.

Step 3: Open HiView and you can view the image, and the photos/videos will be automatically saved in the computer disk.


Note: When connecting the computer to view the image, please make sure to use the line with data transmission function)

4. Work with PC

4.3-USB function introduction.

Select one of the following 3 modes in the USB option of the Settings to enter the mode by default after the microscope is connected to the computer.

- ①**Memory:** Enter the microscope memory from the computer side and view and process the pictures and video files in the current TF card, only the Power button and the Light Wheel are available.
 - ②**PC Camera:**After selecting this mode, you can get the image with larger magnification through the computer camera software. and only the Power button and the Light Wheel are available.
 - ③**Charging only:**When the microscope power is insufficient, it can be charged and used at the same time.
- (***Note:**You can also use your daily phone charger, computer, or power bankto charge the microscope))

5.Specifications

Item	Parameter
Product Name	DM4D
Maximum Video Pixels	1080P
CMOS Size	1/6.5'
Picture Resolution	1920×1080, 1280×720, 640×480
Video Resolution	1080P,720P, VGA
Picture format	JPG
Video Format	AVI
Magnification	500/1000
Lighting	8LED (brightness adjustable)
Focus range	10-40mm (can view distant view)
PC resolution and image transfer rate	640*320/25fps, 1280*720/25fps
PC Operating System	Windows xp,Win7, Win8.1, Win10and above Mac OS x 12 and above
Supported Languages	12 languages
Interface and Signal Transmission Method	USB 2.0
Operating current	420ma
Power	2.1W

6.Warning

- TF card need to choose the transmission speed of C10, if you use a card with a transmission speed lower than C10, the use of pictures, video storage may be unsuccessful and other malfunctions.
- Do not place the microscope in direct sunlight, please avoid rain or humid environment, so as not to cause equipment failure.
- Do not place the microscope near the heat source or fire source, such as electric stove, microwave oven, oven, stove and other environments that may produce high temperature.
- Do not modify the battery, insert foreign objects, immerse in water or other objects to avoid overheating, fire or explosion.
- Please avoid storing the microscope in a place where there is a risk of falling, which may cause damage to the microscope's internal lens deviation or other components, as well as other irreparable damage.