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

Model: HM020

User Manual V3.0


Please read this user manual carefully and keep it for future reference.

Contents

1. Parts Description.....	03
Functions introduction.....	03
2. Interface Introduction.....	05
Menu mode.....	05
Function introduction.....	06
3. Work with PC.....	07
Software download.....	07
Operating steps.....	08
USB Mode Introduction.....	09
4. Specifications.....	10
5. Warm Tips.....	11

1.Parts Description

Notice: The picture is for reference only, please refer to the actual object

• Functions introduction

Part	Function
1	Up button
2	Menu button
3	Down button
4	Screen
5	Power/Photo/Confirm button
6	Indicator light
7	Type-C interface
8	Reset button
9	TF card slot
10	LED lights
11	Camera
12	Transparent cover

Up button: In Settings or albums, it's select button; When previewing the image, it can be digitally zoom in.

Menu button: Short press to enter the Menu interface.

Down button: In Settings or albums, it's select button. In preview mode, short press Down button you can switch between white and purple light.

Screen: 2.0-inch IPS screen.

Power/Photo/Confirm button: Long press for 3 seconds to turn on/off; In Preview mode, press to take photo; In Setting mode, press to confirm.


Indicator light: After powering on, the indicator is long bright blue. After powering off, the indicator is off; The indicator light is red long after charging, green long after full.

Type-C interface: Standard USB2.0 interface, this interface supports charging and data transfer.

Reset button: When the microscope is abnormal, poke this hole with a fine pointed object to force the power off.

TF card slot: Support maximum 32G memory card (Note: please use C10 high-speed card)

LED Lights: There are two sets of lights: white light and UV light, which can be switched.


2.Interface Introduction


(Note: Before using the microscope, to ensure better storage of image files, please insert a TF card into the magnifier)

• Menu mode

Short press the Menu button on the microscope, you can enter the Menu interface, as shown in the figure below:

When viewing files in the album, long-press the menu button to enter the Menu interface, short-press the menu button to Delete Current/ Delete All files, as shown in the figure below:

• Function introduction

Playback: You can view the pictures that have been taken.

Image resolution: 16M, 12M, 10M, 8M, 5M, 3M, 2MHD

Language: Set Microscope language.(English, Portuguese, Simplified Chinese, Japanese, French, German, Korean, Italian, Dutch, Russian, Spanish, Arabic.)

Brightness: Set 8 levels of screen brightness.

Auto Power-off: After selecting auto shutdown time, the microscope will shut down automatically when no one is operating; no auto shutdown, please select "Shut".

Screen savers: Set the screen to extinguish.

Format: This operation will delete all data in the TF card. (Note: It will not be recovered after deletion, please operate carefully)

Default Setting: Restore factory settings.

Version: Check the current firmware version.

3.Work with PC

• Software download

Note: The following 3 methods can be downloaded, please choose the appropriate version according to your computer system.

- ① Please download the applicable software from the following website:
<http://soft.hvscam.com> (Suitable for: Windows 7, Windows 8.1, Windows 10, Windows 11 and Mac))

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
- ② For Windows, search for “HiView” in “Microsoft Store” () ; you can also use Camera () to view. (Suitable for: Windows 8, Windows 10, Windows 11 and above)
- ③ For Mac OS x 10.5 or above, search for “HiView” in “App Store”, you can also use Photo Booth () to view.


• Operating steps

Step 1. Connect the microscope to the computer via a Type-C cable.

Step 2. Long press the microscope power button to turn on, and the computer side of the successful connection screen pop-up window: Memory/PC Camera/Charging Mode three modes, please select the "PC Camera" mode.


Step 3. Run Hiview, you can watch the image, at this time the photo / video files will be automatically saved in the computer disk. (Note: When connecting to a computer to watch images, please make sure you are using a cable with data transfer function)

• USB Mode Introduction

When the device is connected to the computer through a USB cable, the device can be used in Memory /PC Camera / Charging mode.

- ① Memory : Enter the microscope memory from the computer side, you can view and process the current picture files in the TF card, at this time the microscope only Power button and Down button (adjust lights) available.
- ② PC Camera: After selecting this mode, you can get a larger magnification of the image through the computer camera software, at this time the microscope only Power button and Down button (adjust lights) available.
- ③ Charging Mode: When the microscope is low on power, you can use the microscope while charging, use Type-C data cable to connect the microscope and the computer side, the microscope is turned on, select Charging Mode, you can enter the charging state.

(Note: you can also use daily cell phone charger, computer, power bank to charge the microscope)

4.Specifications

Item	Specifications
Product Name	HM020
Screen	2.0-inch IPS screen
CMOS Size	1/4"
Photo Resolution	16M, 12M, 10M, 8M, 5M, 3M, 2MHD
Picture Format	jpg
Illuminant	8 LED lights (four white lights and four purple lights)
Focusing Mode	Manual
Digital Zoom	Supported
White balance	Auto
Exposure	Auto
PC Operating System	Win xp,Win7, Win8.1, Win10, Mac OS x 10.5 or higher
Lens Structure	4P+IR
Aperture	F4.5
Lens angle of view	20°
Interface and Signal Transmission Mode	Type C/usb2.0
Operating Temperature	-10°C ~ +55°C
Operating Humidity	30%~85%Rh
Power Dissipation	260mA,1.3W
Battery capacity	500mAH, 3.7V
Working time	1.5H
Charging time	1H
Language	English, Portuguese, Simplified Chinese, Japanese, French, German, Korean, Italian, Dutch, Russian, Spanish, Arabic.

5. Warm Tips

- TF card with a minimum transfer speed of C10 is recommended. If a card with a transfer speed lower than C10 is used, there may be issues such as unsuccessful storage of photos during use.
- Please avoid placing the microscope in direct sunlight and refrain from exposing it to rain or humid environments to prevent equipment malfunctions.
- Do not place the microscope near a source of heat or fire, such as an electric stove, microwave oven, oven, stove, etc., where high temperatures can be generated.
- Do not modify the battery, insert foreign objects, immerse in water or expose to other substances, in order to prevent overheating, ignition, or explosion.
- Avoid storing the microscope in a place where there is a risk of falling, which may cause lens deviation or damage to other components inside the microscope, as well as other irreparable damage.